

# SELECT SIEVES

## Actions Needed to Specify Sieves

THE FIRST STEP IS TO SELECT THE MEDIA;

- WIRE MESH
- PERFORATED PLATE
- ELECTROFORMED

### 01 SELECT THE FRAME DIAMETER

8 inch - 12 inch

---



### 02 SELECT FRAME HEIGHT


Full or 1/ Height

---

### 03 SELECT FRAME MATERIAL

Brass - SS - Acrylic

---



### 04 BELECT OPENING SIZE

2 to 20 Micron Electro formed -25  
micron to 1MM

---

### 05 SELECT RATING

- Compliance
  - Inspection
  - Calibration
  - Mid-Point
- 

